

This is a special file, named RPTHEAD.TXT, in the directory of a method which allows you to customize the report header page.
It can be used to identify the laboratory which uses the method.

This file is printed on the first page with the report styles:

Header+Short, GLP+Short, GLP+Detail, Short+Spec, Detail+Spec, Full

```

      XXXX  XXX
    XX  XX  XX
  XX      XX      XXXXX  XXX XX
  XX      XX XXX  XX    X  XX X XX
  XX    X  XXX XX  XXXXXXXX  XX X XX
    XX  XX  XX  XX  XX      XX  XX
      XXXX  XXX  XXX  XXXXX  XXX  XXX

```

```

  XXXXXX      X              X      XX
XX    X  XX              XX
XX      XXXXX  XXXXX  XXXXX  XXX      XXXX  XX XXX
  XXXXX  XX      X  XX      XX      XX  XX  XXX XX
      XX  XX      XXXXXX  XX      XX      XX  XX  XX XX
X    XX  XX XX  X  XX      XX XX      XX      XX  XX  XX XX
XXXXXX      XXX  XXXXX X  XXX      XXXX      XXXX  XX  XX

```

```

                                     X
  XX XXX  XXXXX  XX XXX  XXXX  XX XXX  XXXXX
  XXX XX  XX    X  XX  XX  XX  XX  XXX XX  XX
  XX      XXXXXXXX  XX  XX  XX  XX  XX      XX
  XX      XX      XXXXX  XX  XX  XX      XX XX
XXXXX      XXXXX  XX      XXXX  XXXX      XXX
                XXXX

```

```

  XXX              XXX
  XX              XX
  XX      XXXXX  XXXXX  XX      XXXXX  XX XXX
  XX XXX  XX    X      X  XXXXX  XX    X  XXX XX
  XXX XX  XXXXXXXX  XXXXXXXX  XX  XX  XXXXXXXX  XX
  XX  XX  XX      X  XX  XX  XX  XX      XX
  XXX  XXX  XXXXX  XXXXX X  XXXX X  XXXXX  XXXX

```

```

  X              XXX              X
  XX              XX              XX
XXXXX  XXXXX  XXX XX  XX XXX  XX  XXXXX  XXXXX  XXXXX
XX    XX    X  XX X XX  XX  XX  XX      X  XX    XX    X
XX    XXXXXXXX  XX X XX  XX  XX  XX  XXXXXXXX  XX  XXXXXXXX
XX XX  XX      XX  XX  XXXXX  XX  X  XX  XX  XX  XX
  XXX  XXXXX  XXX  XXX  XX      XXXX  XXXXX X  XXX  XXXXX
                XXXX

```

```
=====
Acq. Operator   :                               Seq. Line :    4
Acq. Instrument : Kortemmelab HPLC              Location  : Vial 3
Injection Date  : 11/12/2020 8:06:33 PM          Inj       :    1
                                           Inj Volume : 100.0 µl
Sequence File   : C:\Chem32\1\DATA\20201112_GSP1_F28V_F28Y_S75 2020-11-12 18-26-05\20201112_
                  GSP1_F28V_F28Y_S75.S
Method          : C:\CHEM32\1\DATA\20201112_GSP1_F28V_F28Y_S75 2020-11-12 18-26-05\CJM_S75_
                  RUN_ISOCRATIC.M (Sequence Method)
Last changed    : 11/12/2020 5:53:30 PM
Method Info     : S75 analytical
=====
```

```
=====
Module                                     Type   Firmware rev.   Serial number
-----|-----|-----|-----
1200 Autosampler SL                       G1329B A.06.54 [003] DE64155932
1200 Multiple Wavelength Detector          G1365D B.06.72 [0002] DE64256327
1100/1200 Quaternary Pump                 G1311A A.06.32 [011] DE62971812
1200 Sample Thermostat                    G1330B n/a          DEBAK15882
=====
```

Software Revision: Rev. B.04.03 [16] Copyright © Agilent Technologies

```
=====
                        Column(s)
=====
```

```
Column Description : Sephadex 75, 10/300
Serial#            : 10108795
Product#           : 17-5174-01      Batch# :
Diameter           : 4.6 mm          Length : 150.0 mm
Particle size      : 5.0 µm          Void volume : 60.0 %
# Injections       : 68
Maximum Pressure   : 18.0 bar         Maximum pH : 9.0
Minimum pH         : 2.0
Maximum Temperature: 60.0 °C
Comment           :
=====
```

```
Instrument Conditions :      At Start          At Stop
Pressure             :          0.0             15.7 bar
Flow                 :          0.000           0.800 ml/min
```

Detector Lamp Burn Times: Current On-Time Accumulated On-Time




```
Solvent Description :
PMP1 , Solvent A    : ddH2O
PMP1 , Solvent B    : EtOH
PMP1 , Solvent C    : Tris NaCl
PMP1 , Solvent D    :
```

Run Logbook

12 Nov 20 08:40 PM

Logbook File:C:\Chem32\...Y_S75 2020-11-12 18-26-05\20201112_PE64_F28Y.D\RUN.LOG

Module	# Event Message	Time	Date
Method	Method started: line# 4 vial# 3 inj# 1	20:06:23	11/12/20
Method	Instrument running sample Vial 3	20:06:23	11/12/20
ALS	Air temperature (tray) = 4.0 °C	20:09:45	11/12/20
PUMP	Pressure = 0.0 bar	20:09:45	11/12/20
PUMP	Flow = 0.000 ml/min	20:09:45	11/12/20
PUMP	Pressure = 15.7 bar	20:39:45	11/12/20
Method	Instrument run completed	20:39:48	11/12/20
Method	Saving Method CJM_S75_RUN_ISOCRATIC.M	20:39:49	11/12/20
Method	Saving Method RUN.M	20:39:51	11/12/20
CP Macro	Analyzing rawdata 20201112_PE64_F28Y.D	20:39:51	11/12/20
CP Macro	Can't load MWD1 D, Sig=230,8 Ref=360,100	20:39:51	11/12/20
CP Macro	Signal used in Calib. Table (Signal Details)>	20:39:51	11/12/20



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: MWD1 A, Sig=230,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	12.445	BV	1.3973	4398.21484	41.73980	4.6266
2	15.658	VB	1.0135	9.06645e4	1407.46338	95.3734

Totals : 9.50627e4 1449.20318

Signal 2: MWD1 B, Sig=280,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	12.391	BV	1.3710	625.86176	5.37304	3.1756
2	15.656	VB	0.9765	1.90826e4	306.45151	96.8244

Totals : 1.97085e4 311.82454

Signal 3: MWD1 C, Sig=260,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	12.380	BV	1.4756	564.66089	4.49499	3.6994
2	15.651	VB	0.9803	1.46987e4	234.16295	96.3006

Totals : 1.52634e4 238.65794

=====
*** End of Report ***