

This is a special file, named RPTHEAD.TXT, in the directory of a method which allows you to customize the report header page.

It can be used to identify the laboratory which uses the method.

This file is printed on the first page with the report styles:

Header+Short, GLP+Short, GLP+Detail, Short+Spec, Detail+Spec, Full

```

      XXXX   XXX
    XX  XX   XX
  XX      XX      XXXXX   XXX XX
  XX      XX XXX   XX    X   XX X XX
  XX    X   XXX XX  XXXXXXXX  XX X XX
    XX  XX   XX  XX  XX      XX   XX
      XXXX   XXX   XXX  XXXXX   XXX   XXX

```

```

  XXXXXX      X              X      XX
XX    X   XX              XX
XX      XXXXX   XXXXX   XXXXX   XXX      XXXX   XX XXX
  XXXXX   XX      X   XX      XX      XX  XX   XXX XX
      XX   XX      XXXXXX   XX      XX      XX  XX   XX  XX
X    XX   XX XX  X   XX   XX XX   XX      XX  XX   XX  XX
XXXXXX      XXX   XXXXX X   XXX   XXXX   XXXX   XX  XX

```

```

                                     X
XX XXX   XXXXX   XX XXX   XXXX   XX XXX   XXXXX
  XXX XX  XX    X   XX  XX   XX  XX   XXX XX   XX
  XX      XXXXXXXX  XX  XX   XX  XX   XX      XX
  XX      XX      XXXXX   XX  XX   XX      XX  XX
XXXXX      XXXXX   XX      XXXX   XXXX      XXX
                XXXX

```

```

  XXX              XXX
  XX              XX
  XX      XXXXX   XXXXX   XX      XXXXX   XX XXX
  XX XXX  XX    X      X   XXXXX   XX    X   XXX XX
  XXX XX  XXXXXXXX  XXXXXXX  XX  XX   XXXXXXXX  XX
  XX  XX  XX      X   XX  XX  XX   XX      XX
  XXX   XXX  XXXXX   XXXXX X   XXXX X   XXXXX   XXXX

```

```

  X              XXX              X
  XX            XX              XX
XXXXX   XXXXX   XXX XX   XX XXX   XX   XXXXX   XXXXX   XXXXX
  XX   XX    X  XX X XX   XX  XX   XX      X   XX   XX    X
  XX   XXXXXXXX  XX X XX   XX  XX   XX   XXXXXXXX  XX   XXXXXXXX
  XX XX  XX      XX  XX   XXXXX   XX   X   XX   XX  XX   XX
    XXX   XXXXX   XXX   XXX  XX      XXXX   XXXXX X   XXX   XXXXX
                XXXX

```

```
=====
Acq. Operator   :                               Seq. Line :    3
Acq. Instrument : Kortemmelab HPLC              Location  : Vial 2
Injection Date  : 11/17/2020 11:07:01 PM        Inj       :    1
                                           Inj Volume : 100.0 µl
Sequence File   : C:\Chem32\1\DATA\20201117_GSP1_S75_WITH_STANDARD 2020-11-17 21-00-45\
                  20201117_GSP1_S75_WITH_STANDARD.S
Method          : C:\CHEM32\1\DATA\20201117_GSP1_S75_WITH_STANDARD 2020-11-17 21-00-45\CJM_
                  S75_RUN_ISOCRATIC.M (Sequence Method)
Last changed    : 11/12/2020 5:53:30 PM
Method Info     : S75 analytical
=====
```

```
=====
Module                                     Type   Firmware rev.   Serial number
-----|-----|-----|-----
1200 Autosampler SL                       G1329B A.06.54 [003] DE64155932
1200 Multiple Wavelength Detector          G1365D B.06.72 [0002] DE64256327
1100/1200 Quaternary Pump                 G1311A A.06.32 [011] DE62971812
1200 Sample Thermostat                   G1330B n/a          DEBAK15882
=====
```

Software Revision: Rev. B.04.03 [16] Copyright © Agilent Technologies

```
=====
                        Column(s)
=====
```

```
Column Description : Sephadex 75, 10/300
Serial#            : 10108795
Product#           : 17-5174-01      Batch# :
Diameter           : 4.6 mm          Length : 150.0 mm
Particle size      : 5.0 µm          Void volume : 60.0 %
# Injections       : 73
Maximum Pressure   : 18.0 bar         Maximum pH : 9.0
Minimum pH         : 2.0
Maximum Temperature: 60.0 °C
Comment           :
=====
```

```
Instrument Conditions :      At Start          At Stop
Pressure              :          0.0             15.7 bar
Flow                  :          0.000           0.800 ml/min
```

Detector Lamp Burn Times: Current On-Time Accumulated On-Time


```
Solvent Description :
PMP1 , Solvent A    : ddH2O
PMP1 , Solvent B    : EtOH
PMP1 , Solvent C    : Tris NaCl
PMP1 , Solvent D    :
```

Run Logbook

17 Nov 20 11:40 PM

Logbook File:C:\Chem32\...ANDARD 2020-11-17 21-00-45\20201117_STANDARD.D\RUN.LOG

Module	# Event Message	Time	Date
Method	Method started: line# 3 vial# 2 inj# 1	23:06:50	11/17/20
Method	Instrument running sample Vial 2	23:06:51	11/17/20
ALS	Air temperature (tray) = 4.0 °C	23:10:12	11/17/20
PUMP	Pressure = 0.0 bar	23:10:12	11/17/20
PUMP	Flow = 0.000 ml/min	23:10:12	11/17/20
PUMP	Pressure = 15.7 bar	23:40:13	11/17/20
Method	Instrument run completed	23:40:16	11/17/20
Method	Saving Method CJM_S75_RUN_ISOCRATIC.M	23:40:17	11/17/20
Method	Saving Method RUN.M	23:40:18	11/17/20
CP Macro	Analyzing rawdata 20201117_STANDARD.D	23:40:18	11/17/20
CP Macro	Can't load MWD1 D, Sig=230,8 Ref=360,100	23:40:18	11/17/20
CP Macro	Signal used in Calib. Table (Signal Details)>	23:40:18	11/17/20



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: MWD1 A, Sig=230,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	11.194	BV	0.6411	1.66567e5	3169.65942	19.7937
2	11.910	VV	0.9188	2.20131e5	2851.11279	26.1589
3	14.414	VV	1.0865	2.09873e5	2636.93359	24.9400
4	16.911	VV	0.3826	5.42382e4	1698.14209	6.4453
5	16.986	VV	0.0784	9403.99609	1703.74353	1.1175
6	17.066	VB	0.4234	5.94183e4	1701.44397	7.0609
7	23.382	BV	0.6479	1.20516e5	2218.11084	14.3213
8	26.204	VBA	0.9797	1366.15295	19.20987	0.1623

Totals : 8.41513e5 1.59984e4

Signal 2: MWD1 B, Sig=280,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	11.194	BV	0.5170	7.63285e4	2187.92261	23.9841
2	11.950	VV	0.8775	9.01830e4	1474.59827	28.3375
3	14.417	VV	1.0986	6.56622e4	843.13531	20.6325
4	17.025	VV	0.1270	49.63631	5.21702	0.0156
5	18.773	VV	1.7583	2.26215e4	160.11418	7.1082
6	23.378	VB	1.4307	5.59603e4	511.37949	17.5840
7	26.233	BV	1.4035	7440.69189	66.49407	2.3380

Totals : 3.18246e5 5248.86095

Signal 3: MWD1 C, Sig=260,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	11.187	BV	0.5610	6.26845e4	1614.30994	20.8747
2	11.926	VV	0.9398	6.03231e4	897.91577	20.0883
3	14.418	VV	1.5137	7.45617e4	657.14923	24.8299
4	18.824	VV	2.0473	4.57509e4	334.57715	15.2356
5	23.378	BB	0.9620	3.63434e4	523.47638	12.1028
6	26.300	BV	1.9355	1.91349e4	119.15430	6.3721
7	29.551	VBA	0.9699	1491.28284	18.20095	0.4966

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
----- ----- ----- ----- ----- ----- -----						
Totals :				3.00290e5	4164.78372	

=====
*** End of Report ***