

This is a special file, named RPTHEAD.TXT, in the directory of a method which allows you to customize the report header page.
It can be used to identify the laboratory which uses the method.

This file is printed on the first page with the report styles:

Header+Short, GLP+Short, GLP+Detail, Short+Spec, Detail+Spec, Full

```

      XXXX   XXX
    XX  XX   XX
  XX      XX      XXXXX   XXX XX
  XX      XX XXX   XX    X   XX X XX
  XX      X   XXX XX XXXXXXXX XX X XX
    XX XX   XX  XX  XX      XX   XX
      XXXX   XXX   XXX   XXXXX   XXX   XXX

```

```

  XXXXXX      X              X      XX
XX   X   XX              XX
XX      XXXXX   XXXXX   XXXXX   XXX   XXXX   XX XXX
  XXXXX   XX      X   XX   XX   XX   XX   XX   XXX XX
      XX   XX   XXXXXX   XX   XX   XX   XX   XX   XX
X   XX   XX XX   X   XX   XX XX   XX   XX   XX   XX
XXXXXX      XXX   XXXXX X   XXX   XXXX   XXXX   XX   XX

```

```

                                     X
XX XXX   XXXXX   XX XXX   XXXX   XX XXX   XXXXX
  XXX XX  XX    X   XX  XX   XX  XX   XXX XX   XX
  XX      XXXXXXXX   XX  XX   XX  XX   XX      XX
  XX      XX      XXXXX   XX  XX   XX      XX XX
XXXXX      XXXXX   XX      XXXX   XXXX      XXX
                XXXX

```

```

  XXX              XXX
  XX              XX
  XX      XXXXX   XXXXX   XX   XXXXX   XX XXX
  XX XXX  XX    X      X   XXXXX   XX    X   XXX XX
  XXX XX  XXXXXXXX   XXXXXXX   XX  XX   XXXXXXXX   XX
  XX  XX  XX      X   XX  XX  XX   XX      XX
  XXX   XXX   XXXXX   XXXXX X   XXXX X   XXXXX   XXXX

```

```

  X              XXX              X
  XX              XX              XX
XXXXXX   XXXXX   XXX XX   XX XXX   XX   XXXXX   XXXXX   XXXXX
  XX   XX    X   XX X XX   XX  XX   XX      X   XX   XX    X
  XX   XXXXXXXX   XX X XX   XX  XX   XX   XXXXXXXX   XX   XXXXXXXX
  XX XX  XX      XX  XX   XXXXX   XX   X   XX   XX XX   XX
    XXX   XXXXX   XXX   XXX   XX      XXXX   XXXXX X   XXX   XXXXX
                XXXX

```

```
=====
Acq. Operator   :                               Seq. Line :    7
Acq. Instrument : Kortemmelab HPLC              Location  : Vial 4
Injection Date  : 11/18/2020 1:21:54 AM          Inj       :    1
                                           Inj Volume : 100.0 µl
Sequence File   : C:\Chem32\1\DATA\20201117_GSP1_S75_WITH_STANDARD 2020-11-17 21-00-45\
                  20201117_GSP1_S75_WITH_STANDARD.S
Method          : C:\CHEM32\1\DATA\20201117_GSP1_S75_WITH_STANDARD 2020-11-17 21-00-45\CJM_
                  S75_RUN_ISOCRATIC.M (Sequence Method)
Last changed    : 11/12/2020 5:53:30 PM
Method Info     : S75 analytical
=====
```

```
=====
Module                                Type   Firmware rev.   Serial number
-----|-----|-----|-----
1200 Autosampler SL                  G1329B A.06.54 [003] DE64155932
1200 Multiple Wavelength Detector    G1365D B.06.72 [0002] DE64256327
1100/1200 Quaternary Pump            G1311A A.06.32 [011] DE62971812
1200 Sample Thermostat               G1330B n/a          DEBAK15882
=====
```

Software Revision: Rev. B.04.03 [16] Copyright © Agilent Technologies

=====

Column(s)

=====

```
Column Description : Sephadex 75, 10/300
Serial#            : 10108795
Product#           : 17-5174-01      Batch# :
Diameter           : 4.6 mm          Length : 150.0 mm
Particle size      : 5.0 µm          Void volume : 60.0 %
# Injections       : 77
Maximum Pressure   : 18.0 bar         Maximum pH : 9.0
Minimum pH         : 2.0
Maximum Temperature: 60.0 °C
Comment           :
=====
```

```
Instrument Conditions :      At Start          At Stop
Pressure             :          0.0              15.6 bar
Flow                 :          0.000            0.800 ml/min
```

Detector Lamp Burn Times: Current On-Time Accumulated On-Time




```
Solvent Description :
PMP1 , Solvent A    : ddH2O
PMP1 , Solvent B    : EtOH
PMP1 , Solvent C    : Tris NaCl
PMP1 , Solvent D    :
```

Run Logbook

18 Nov 20 01:55 AM

Logbook File:C:\Chem32\...H_STANDARD 2020-11-17 21-00-45\20201117_F28V.D\RUN.LOG

Module	# Event Message	Time	Date
Method	Method started: line# 7 vial# 4 inj# 1	01:21:44	11/18/20
Method	Instrument running sample Vial 4	01:21:44	11/18/20
ALS	Air temperature (tray) = 4.0 °C	01:25:06	11/18/20
PUMP	Pressure = 0.0 bar	01:25:06	11/18/20
PUMP	Flow = 0.000 ml/min	01:25:06	11/18/20
PUMP	Pressure = 15.6 bar	01:55:07	11/18/20
Method	Instrument run completed	01:55:09	11/18/20
Method	Saving Method CJM_S75_RUN_ISOCRATIC.M	01:55:11	11/18/20
Method	Saving Method RUN.M	01:55:14	11/18/20
CP Macro	Analyzing rawdata 20201117_F28V.D	01:55:14	11/18/20
CP Macro	Can't load MWD1 D, Sig=230,8 Ref=360,100	01:55:14	11/18/20
CP Macro	Signal used in Calib. Table (Signal Details)>	01:55:14	11/18/20



=====
Area Percent Report
=====

Sorted By : Signal
Multiplier: : 1.0000
Dilution: : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: MWD1 A, Sig=230,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	11.158	BV	0.5679	3211.05884	77.73979	8.7216
2	12.558	VV	1.3550	7013.91992	69.40488	19.0507
3	15.251	VB	0.8844	1.74260e4	287.80521	47.3312
4	22.167	BV	0.7096	8853.43262	183.26131	24.0470
5	27.618	BB	0.4919	312.74966	7.76500	0.8495

Totals : 3.68172e4 625.97618

Signal 2: MWD1 B, Sig=280,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	0.110	BB	0.0368	8.55713	4.01798	0.1717
2	11.201	BV	0.6316	670.01050	14.27353	13.4459
3	12.016	VV	1.1982	960.59125	9.52812	19.2774
4	15.251	VB	0.8491	3343.83643	58.12901	67.1049

Totals : 4982.99530 85.94864

Signal 3: MWD1 C, Sig=260,8 Ref=360,100

Peak #	RetTime [min]	Type	Width [min]	Area [mAU*s]	Height [mAU]	Area %
1	0.113	BB	0.0226	1.77670	1.95279	0.0415
2	11.226	BB	0.9465	1257.42615	17.28959	29.3542
3	15.251	BB	0.8694	2657.47461	44.58674	62.0378
4	22.185	BV	1.1028	366.95786	4.26911	8.5665

Totals : 4283.63531 68.09823

=====
*** End of Report ***